

## <sup>13</sup>C Breast Array

### MR Coils Made-to-Measure

- 8 channel array:
  - 2 channel Tx/Rx
  - 6 receive elements
- supports parallel imaging
- use of <sup>1</sup>H body coil allowed
- CE certified as medical device
- UL 94 approved
- enhanced patient comfort by ergonomically shaped housing and cushion support set
- easy handling due to low coil weight



### Specifications

approvals	Medical Device, CE according to Council Directive 93/42/EEC
$B_0$ -field strength	3 T
housing dimensions	37 x 54 x 17.5 mm
weight	ca. 4.7 kg (coil only)