

^{13}C Transmit Coil*

MR Coils Made-to-Measure

- ^{13}C transmit coil
- transmit-only (Tx-only), linear polarization
- Helmholtz-like element arrangement, 2 loops (300 x 300 mm²)
- optimized for prostate, abdominal excitation with reduced efficiency
- detachable top for easy patient positioning
- maximum transmit power: 6.5 kW (MR system specific)
- interfacing to MR system
- use of ^1H body coil allowed



Please contact us for availability on your MR system.

Specifications

* Regulatory requirements for medical products will vary by country and MR system. Please contact us at info@rapidbiomed.de or info@rapidmri.com (USA, Canada, and South America) to determine approval status for products mentioned on this product sheet.

B_0 -field strength	3 T
housing dimensions	45 cm x 53 cm x 41 cm (length x width x height)
weight	ca. 7 kg