

¹³C Breast Array

MR Coils Made-to-Measure

- 8 channel array:
 - 2 channel Tx/Rx
 - 6 receive elements
- supports parallel imaging
- use of ¹H body coil allowed
- CE certified as medical device
- UL 94 approved
- enhanced patient comfort by ergonomically shaped housing and cushion support set
- easy handling due to low coil weight



Specifications

* Regulatory requirements for medical products will vary by country and MR system. Please contact us at info@rapidbiomed.de or info@rapidmri.com (USA, Canada, and South America) to determine approval status for products mentioned on this product sheet.

B_0 -field strength	3 T
housing dimensions	37 x 54 x 17.5 mm
weight	ca. 4.7 kg (coil only)