

## 7 T X Nucleus Cardiac Coils \*

### MR Coils Made-to-Measure



- dual tuned surface coil for cardiac studies
- rigid housing
- x-nucleus channel with linear polarization
- proton channel with quadrature polarization
- available in any combination of nuclei required for your studies (i.e.  $^{31}\text{P}/^1\text{H}$ ,  $^{23}\text{Na}/^1\text{H}$ ,  $^{13}\text{C}/^1\text{H}$ , ...).
- for prone or supine patient positioning
- anatomically shaped coil housing
- individually adaptable to most established 7 T human MR systems

## Specifications

\* Regulatory requirements for medical products will vary by country and MR system. Please contact us at [info@rapidbiomed.de](mailto:info@rapidbiomed.de) or [info@rapidmri.com](mailto:info@rapidmri.com) (USA, Canada, and South America) to determine approval status for products mentioned on this product sheet.

$B_0$ -field strength	7 T
dimensions	49 x 30 x 5.5 cm
weight	2.9 kg