


^{13}C Paddle Coil *

MR Coils Made-to-Measure



^{13}C transmit coil

- ^{13}C multi-purpose Rx array
- two coil paddles
- 8 x ^{13}C Rx channels (4 channels each)
- flexible coil positioning
- use of ^1H body coil allowed
- dedicated ^{13}C transmit coil required

Individually adaptable to most established MR systems with minimum 8 receive channels. Please contact us for availability on your MR system.

Specifications

* Regulatory requirements for medical products will vary by country and MR system. Please contact us at info@rapidbiomed.de or info@rapidmri.com (USA, Canada, and South America) to determine approval status for products mentioned on this product sheet.

B_0 -field strength	3 T
dimensions of coil housing	24 x 21 x 6 cm ³ each paddle
weight	2.6 kg