

^{13}C Abdominal Arrays

MR Coils Made-to-Measure



Specifications

- ^{13}C receive-only
- up to 16 independent receive channels
- supports parallel imaging
- dedicated ^{13}C transmit coil required
- use of ^1H body coil allowed
- B_0 -field strength 3 T
- flexible housing option
- individually adaptable to different MR systems

Regulatory requirements for medical devices will vary by country and MR system.

Please, contact us at info@rapidbiomed.de or info@rapidmri.com (USA, Canada, and South America) to determine approval status for products mentioned on this product sheet.

Pre-clinical results from pig studies




Image Courtesy: Christoffer Laustsen, Aarhus University, Denmark