


## $^{13}\text{C}$ Breast Array

MR Coils Made-to-Measure



### Specifications

- phased array for  $^{13}\text{C}$  breast imaging
- 8 channel array:
  - 2 channel Tx/Rx
  - 6 receive elements
- supports parallel imaging
- use of  $^1\text{H}$  body coil allowed
- $B_0$ -field strength: 3 T
- ergonomically shaped coil housing and cushions
- low profile housing in favour of patient positioning
- approx. housing dimensions: 37 x 54 x 20 cm
- approx. weight: 4.7 kg (coil only)
- individually adaptable to different MR systems

Regulatory requirements for medical devices will vary by country and MR system.

Please, contact us at [info@rapidbiomed.de](mailto:info@rapidbiomed.de) or [info@rapidmri.com](mailto:info@rapidmri.com) (USA, Canada, and South America) to determine approval status for products mentioned on this product sheet.