


3 T Dual Tuned Head Coils

MR Coils Made-to-Measure



Data Courtesy: Jan Ruff, Siemens AG

Specifications

- volume coils for studies of the human head
- both channels equipped with quadrature polarization
- capable of performing ^1H decoupled experiments when applicable
- available for ^{31}P , ^{23}Na and ^{13}C combined with ^1H (other non-proton nuclei on demand)
- B_0 -field strength 3 T
- housing dimensions: inner diameter 26.5 cm, housing length 43 cm
- approx. weight: 10 kg
- a pediatric version with adapted dimensions is also available
- individually adaptable to different MR systems

Regulatory requirements for medical devices will vary by country and MR system.

Please, contact us at info@rapidbiomed.de or info@rapidmri.com (USA, Canada, and South America) to determine approval status for products mentioned on this product sheet.