


## $^{13}\text{C}$ Paddle Coil

MR Coils Made-to-Measure



### Specifications

- $^{13}\text{C}$  multi-purpose Rx array
- two coil paddles
- 8 Rx channels (4 channels each)
- flexible coil positioning
- dedicated  $^{13}\text{C}$  transmit coil required
- use of  $^1\text{H}$  body coil allowed
- $B_0$ -field strength 3 T
- approx. housing dimensions: 24 x 21 x 6 cm<sup>3</sup> each paddle
- approx. weight: 2.6 kg
- individually adaptable to different MR systems with minimum 8 receive channels

Regulatory requirements for medical devices will vary by country and MR system.

Please, contact us at [info@rapidbiomed.de](mailto:info@rapidbiomed.de) or [info@rapidmri.com](mailto:info@rapidmri.com) (USA, Canada, and South America) to determine approval status for products mentioned on this product sheet.