


Non-Proton Torso Coils

MR Coils Made-to-Measure



Specification

- non-proton volume coils for human torso
- open design for elbow-room
- housing adapted to magnet bore and patient table
- upper part detachable for patient access
- 4 rung quadrature birdcage
- for B_0 fields up to 7 T
- adaptable to GE, Philips or Siemens MR systems

Add a receive array for SNR performance.

Please, contact us at info@rapidbiomed.de or info@rapidmri.com (USA, Canada, and South America) to determine approval status for products mentioned on this product sheet.