

## 4 Ch Obese Rat Thorax Array

### MR Coils Made-to-Measure




Image Courtesy Detlef Stiller, Boehringer Ingelheim Pharma GmbH & Co. KG Div. Research, Germany

- 4 channels for optimised image quality
- whole thorax coverage by volumetric array
- split coil housing supports
  - easy animal handling
  - reproducible animal positioning
  - fast animal change for maximum throughput
- integrated preamplifiers for optimum sensitivity
- large diameter for obese rats

## Specifications

**No Medical Device! Caution** – The use of the devices described above is limited to investigational use on laboratory animals or other tests that do not involve human subjects.

coil	rat thorax array coil
features	volumetric coil array 4 channels split housing for easy animal handling
$B_0$ field strength	1.5 T – 9.4 T
required equipment	actively decoupled volume transmit resonator with minimum inner diameter 86 mm