

PRODUCT DATASHEET

Dual Tuned Extremity Coil 7 T

Specifications

- dual tuned coil for calf imaging
- Transmit / receive coil
- both channels equipped with quadrature polarization
- available in various combinations of nuclei (e.g. ^{31}P / ^1H , ^{23}Na / ^1H , ^{13}C / ^1H ...)
- approx. housing dimensions of volume resonator: ID 22 cm
- lower leg holder for coil changes without repositioning
- for 3 T and 7 T MR systems
- no medical device
- optional: dedicated phantom with different compartments for reference concentrations



lower leg support
including 5 compartments for
phantom fluids and heel cushion

Regulatory requirements for medical devices will vary by country and MR system.

Please, contact us at info@rapidbiomed.de or info@rapidmri.com (USA, Canada, and South America) to determine approval status for products mentioned on this product sheet.