

PRODUCT DATASHEET

Dual Tuned Cardiac Coil 7 T

Specifications

- dual tuned surface coil for cardiac studies
- X nucleus channel with linear polarization
- proton channel with quadrature polarization
- available in various combinations of nuclei (e.g. ^{31}P / ^1H , ^{23}Na / ^1H , ^{13}C / ^1H , other nuclei on demand)
- approx. housing dimensions: 30 x 49 cm (w x l)
- rigid, curved coil housing
- for prone or supine patient positioning
- no medical device



Regulatory requirements for medical devices will vary by country and MR system.

Please, contact us at info@rapidbiomed.de or info@rapidmri.com (USA, Canada, and South America) to determine approval status for products mentioned on this product sheet.