


PRODUCT DATASHEET

X Nucleus Cardiac Coils

Specifications

- flexible or rigid housings
- available in linear or quadrature polarization or as multi channel receive array
- single tuned or dual tuned
- for prone or supine patient positioning, depending on coil design
- adjustable housings for e.g. pediatric studies
- for 3 T and 7 T MR systems
- nuclei: ^1H , ^{31}P , ^{23}Na , ^{13}C , ^{19}F and other



Regulatory requirements for medical devices will vary by country and MR system.

Please, contact us at info@rapidbiomed.de or info@rapidmri.com (USA, Canada, and South America) to determine approval status for products mentioned on this product sheet.