


PRODUCT DATASHEET

IMCL Customized Setup for Mouse Leg Muscle

Specifications

- dedicated holder system
- anesthesia mask with adjustable bite bar
- anatomically shaped mouse bed
- clamping mechanism for reproduceable positioning of mouse leg muscle and coil

No Medical Device! Caution – The use of the devices described above is limited to investigational use on laboratory animals or other tests that do not involve human subjects.



No Medical Device! Caution – The use of the devices described above is limited to investigational use on laboratory animals or other tests that do not involve human subjects.