

# PRODUCT DATASHEET

## Mouse Holders

### Mouse Holder LHMXGS-00562

- integrated gas anesthesia and stereotactic holder
- outer diameter 71 mm
- length 87 cm

including heart and  
head array



This holder is compatible with the RAPID standard volume transmit coils with an inner diameter of 72 mm and our mouse brain and heart arrays.

### Mouse Holder for Small Volume Coils LHMXG-01017

- integrated gas anesthesia
- outer diameter 34 mm
- length 300 mm

This holder is compatible with the RAPID standard volume mouse coils.

We also offer alternative outer diameters.



**No Medical Device! Caution** – The use of the devices described above is limited to investigational use on laboratory animals or other tests that do not involve human subjects.