

PRODUCT DATASHEET

Rat Holders

Rat Holder LHRXGS-00563

- integrated gas anesthesia and stereotactic holder
- outer diameter 71 mm
- length 87 cm



This holder is compatible with the RAPID standard volume transmit coils with an inner diameter of 72 mm and our rat brain and heart arrays.

Rat Holder for Obese Rats LHRXGS-01107

- integrated gas anesthesia and stereotactic holder
- outer diameter 89 mm
- length 70 cm



This holder is compatible with the RAPID volume coils for obese rats with an inner diameter of 90 mm.


No Medical Device! Caution – The use of the devices described above is limited to investigational use on laboratory animals or other tests that do not involve human subjects.

Rat Holders

Rat Holder for Cardiac Studies LHRXGS-00787

- integrated gas anesthesia and stereotactic holder
- outer diameter 71 mm, length 76 cm
- cardiac coil slides under the head holder for flexible coil positioning

This holder is compatible with the RAPID rat cardiac arrays.



No Medical Device! Caution – The use of the devices described above is limited to investigational use on laboratory animals or other tests that do not involve human subjects.